Jamshidi Evolve™

- en** Jamshidi Evolve™ Bone Marrow Biopsy/Aspiration Needle
- fr** Aiguille d'aspiration/de biopsie Jamshidi Evolve™ pour moelle osseuse
- de** Jamshidi Evolve™ Knochenmarkbiopsie-/Aspirationsnadel
- it** Ago per aspirazione/biopsia del midollo osseo Jamshidi Evolve™
- es** Aguja para aspiración/biopsia de médula ósea Jamshidi Evolve™
- pt** Agulha de aspiração/biópsia de medula óssea Jamshidi Evolve™
- nl** Jamshidi Evolve™-beenmergbiopsie-/beenmergaspiratiennaald
- sv** Jamshidi Evolve™ benmärgsbiopsi-/aspirationsnål
- da** Jamshidi Evolve™ knoglemarvsbiopsi-/aspirationsnål
- fi** Jamshidi Evolve™-luuydinbiopsia-/aspiraationeula
- no** Jamshidi Evolve™-nål til benmargsbiopsi / aspirasjon
- el** Βελόνα βιοψίας/αναρρόφησης μυελού των οστών Jamshidi Evolve™
- ja** Jamshidi Evolve™ 骨髓穿刺 / 吸引針
- zh** Jamshidi Evolve™ 骨髓活检/抽吸針

BD, the BD Logo, and Jamshidi Evolve are trademarks of Becton, Dickinson and Company or its affiliates. © 2020 BD. All rights reserved.

Manufacturer:

Bard Peripheral Vascular, Inc.

1625 West 3rd Street
Tempe, AZ 85281

USA

TEL: 1-480-894-9515

1-800-321-4254

FAX: 1-480-966-7062

1-800-440-5376

www.crbard.com/biopsy

**Authorized Representative
in the European Community:**
BD Switzerland Sàrl
Terre Bonne Park - A4
Route de Crassier 17
1262 Eysins - Switzerland

361-40301

08/2019

STERILE EO

Sterilized Using
Ethylene Oxide

For Single
Use Only

Do Not
Resterilize

Do not use if the product sterile
barrier system or its packaging
is compromised



R_x Only

CAUTION: Federal (USA) law
restricts this device to sale by
or on the order of a physician


BD


Figure, Abbildung, Figura, Afbeelding, Figur, Kuva, Εικόνα, 图, 圖

(4)


(1)


(5)

(2)


(6)


(3)


Figure, Abbildung, Figura, Afbeelding, Figur, Kuva, Εικόνα, 图, 图


(7)


(10)

(8)

(11)

(9)

This instructions for use is designed to serve only as a general guideline. It is not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

Indications for Use

Indicated for use in aspirating bone marrow and for use in obtaining core biopsy samples of bone and/or bone marrow.

Contents

- A – Cap
- B – Stylet
- C – Handle
- D – Cannula
- E – Probe
- F – Probe Guide
- G – Specimen Cradle (included in EJM only)
- H – Depth Stop

Warnings

- Only clinicians who have been properly trained in this procedure, including patient preparation, sample preservation and, if required, the use of image guidance for percutaneous biopsy, should use this device. The instructions are NOT intended to define or suggest surgical protocol. It is the physician's responsibility to determine the appropriate protocol and techniques required based upon the needs of the patient.
- Good medical judgment should be exercised in the proper selection of needle dimensions and use of image guidance based on biopsy location and patient anatomy.
- Good medical judgment should be exercised in considering biopsy on patients who are receiving anticoagulant therapy or who have bleeding disorders.
- Do not use for a sternal procedure. Due to the needle length, internal thoracic organs or blood vessels may be punctured or otherwise damaged.

Cautions

- For single use only. Re-use may result in a non-functional product or contribute to cross contamination.
- Do not re-sterilize.
- Do not use if package is damaged.
- Prior to needle insertion, inspect the needle for a damaged point, bent shaft or other imperfections.
- Applying too much axial force may bend the needle. Angularity beyond 20° may cause the needle to bend or break.
- After use, this product may be a biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state, and federal laws and regulations.

Potential Complications

Potential complications associated with aspiration and core biopsy procedures are site specific and include, but are not limited to: cell seeding, infection, pain, bleeding, hematoma, non-target tissue, adjacent tissue injury, and organ or vessel perforation.

Instructions for Use

Procedure Preparation

1. Determine the desired anatomical site and select the proper needle size and length.

Optional: If using the depth stop, place the depth stop to the desired location on the cannula.

2. Place the patient in an appropriate position for the biopsy procedure. Locate and mark the appropriate needle insertion site. **(Figure 1)**
3. Using sterile technique, prepare the site.
4. Infiltrate the skin, especially the periosteum, with local anesthetic.
5. A skin puncture may be created prior to inserting the needle.

Needle Insertion

6. Hold the handle securely in the palm of your hand and introduce the needle through the skin. **(Figure 2)**
7. Using gentle but firm pressure, slowly advance the needle while rotating the needle in an alternating clockwise-counterclockwise motion. Entrance into the marrow cavity is generally detected by decreased resistance. **(Figure 3)**
8. Rotate the stylet 90° counterclockwise and remove the stylet.

Aspiration

When using the cannula for biopsy, the aspiration should be taken before the biopsy procedure (EJC). When using the cradle for biopsy, the aspiration may be taken before or after the biopsy procedure (EJM).

9. Prepare the syringe per hospital or facility protocol.
10. Attach a syringe to the cannula hub. **(Figure 4)**
11. Apply suction by quickly withdrawing the syringe plunger.
12. Disconnect the syringe and remove the aspirated specimen.
13. If the procedure is complete, slowly remove the cannula using a clockwise-counterclockwise motion and properly dispose of the stylet and cannula.
14. If a biopsy sample will also be collected, use the same location or choose a different location to enter the bone. If a different location is used, follow the instructions for Needle Insertion above. Follow the instructions for Biopsy using Cannula or Biopsy using Cradle below.

Biopsy using Cannula (EJC)

15. After removal of the stylet, the cap may be placed onto the handle to maintain handle comfort during the procedure.
16. Centimeter depth markings and the depth stop may be used for guiding insertion depth and approximating biopsy sample length.
17. With the stylet removed, slowly advance the cannula with a slight clockwise-counterclockwise motion until adequate marrow is obtained.

Note: Sample architecture may be damaged with excessive clockwise-counterclockwise movement.

18. With an angular motion, completely rotate the cannula while applying light downward pressure to allow the tapered tip to help sever the biopsy specimen. **(Figure 5)**

Note: This is a suggested method for sample severing; where applicable, additional rotations may be performed or follow your hospital or facility protocol.

Optional: Attach a syringe to the handle and draw negative pressure. This may assist biopsy sample retrieval.

19. Slowly remove the cannula using a clockwise-counterclockwise motion.

20. To eject the biopsy sample, remove the cap or syringe from the cannula hub as needed. Place the probe guide onto the distal cannula tip. Remove the biopsy sample by inserting the probe into the cannula tip. **(Figure 6)**
21. To collect an additional biopsy sample, verify the cannula is not compromised from its original state. Repeat the applicable steps above.
22. When the procedure is complete, properly dispose of the stylet and cannula.

Biopsy using Cradle (EJM)

15. After removal of the stylet, the cap may be placed onto the handle to maintain handle comfort during the procedure.
16. Centimeter depth markings and the depth stop may be used for guiding insertion depth and approximating biopsy sample length.
17. With the stylet removed, slowly advance the cannula with a slight clockwise-counterclockwise motion until adequate marrow is obtained.

Note: Sample architecture may be damaged with excessive clockwise-counterclockwise movement.

18. Insert the cradle into the cannula.

Optional: Advance only far enough to rest on top of the biopsy sample. Approximate biopsy sample length can be identified at this time. If needed, advance the cannula further inside the bone marrow cavity.

19. Advance the cradle down to the handle base or until resistance is felt. **(Figure 7)**
20. Rotate the cannula handle and cradle together. **(Figure 8)**
21. Remove the cradle from the cannula. **(Figure 9)**

Note: The arrow on the white hub of the cradle indicates the direction of the opening where the biopsy sample will be located. **(Figure 10)**

22. Advance the probe through the cradle hub to remove the biopsy sample. **(Figure 11)**
23. If the procedure is complete, slowly remove the cannula using a clockwise-counterclockwise motion.
24. To collect an additional biopsy sample, verify the cradle is not compromised from its original state and repeat the applicable steps above.
25. To collect an aspiration after biopsy, use the same location or choose a different location to enter the bone. Follow the instructions for Aspiration above.
26. When the procedure is complete, properly dispose of the stylet and cannula.

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